

## Anti-MYOC antibody {FITC} (STJ501851)

| GENERAL INFORMATION    |                     |                                                                                                                                                                      |
|------------------------|---------------------|----------------------------------------------------------------------------------------------------------------------------------------------------------------------|
|                        | Product Type        | Primary antibodies                                                                                                                                                   |
|                        | Applications        | CM/ELISA/ICC/IF/IHC/IP/WB                                                                                                                                            |
|                        | Host / Source       | Rabbit                                                                                                                                                               |
|                        | Reactivity          | Human/Monkey/Mouse/Rat                                                                                                                                               |
| PRODUCT PROPERTIES     |                     |                                                                                                                                                                      |
|                        | Clonality           | Polyclonal                                                                                                                                                           |
|                        | Concentration       | 0.75 µg/µl                                                                                                                                                           |
|                        | Conjugation         | FITC                                                                                                                                                                 |
|                        | Purification        | Affinity Purified                                                                                                                                                    |
|                        | Dilution Range      | WB: 1:500<br>DB: 1:10, 000<br>ELISA: 1:10, 000<br>IP: 1:200<br>IHC: 1:200<br>ICC: 1:200<br>IF: 1:200<br>CFM: 1:200                                                   |
|                        | Formulation         | Contains Tris, HCl/Glycine buffer pH 7.4-7.8, 30% Glycerol and 0.5% BSA, along with cryo-protective agents, Hepes, and long-term preservatives (0.02% Sodium Azide). |
|                        | Isotype             | IgG                                                                                                                                                                  |
|                        | Storage Instruction | Store at-20°C for long term storage. Avoid freeze-thaw cycles.                                                                                                       |
| TARGET INFORMATION     |                     |                                                                                                                                                                      |
|                        | Gene ID             | <a href="#">4653</a>                                                                                                                                                 |
|                        | Gene Symbol         | <a href="#">MYOC</a>                                                                                                                                                 |
|                        | UniProt ID          | <a href="#">MYOC_HUMAN</a>                                                                                                                                           |
| ADDITIONAL INFORMATION |                     |                                                                                                                                                                      |
|                        | Note                | STRICTLY FOR FURTHER SCIENTIFIC RESEARCH USE ONLY (RUO). MUST NOT TO BE USED IN DIAGNOSTIC OR THERAPEUTIC APPLICATIONS.                                              |