

## Anti-KIRREL antibody (200-250) (STJ501885)

### GENERAL INFORMATION

|               |                       |
|---------------|-----------------------|
| Product Type  | Primary antibodies    |
| Applications  | CM/ELISA/IF/IHC/IP/WB |
| Host / Source | Rabbit                |
| Reactivity    | Mouse/Rat             |

### PRODUCT PROPERTIES


|                     |                                                                                                                                                                      |
|---------------------|----------------------------------------------------------------------------------------------------------------------------------------------------------------------|
| Clonality           | Polyclonal                                                                                                                                                           |
| Concentration       | 0.45-0.65 µg/µl                                                                                                                                                      |
| Conjugation         | Unconjugated                                                                                                                                                         |
| Purification        | Affinity Purified                                                                                                                                                    |
| Dilution Range      | WB: 1:250<br>DB: 1:10, 000<br>ELISA: 1:10, 000<br>IP: 1:250<br>IHC: 1:250<br>ICC: 1:250<br>IF: 1:250<br>CFM: 1:250                                                   |
| Formulation         | Contains Tris, HCl/Glycine buffer pH 7.4-7.8, 30% Glycerol and 0.5% BSA, along with cryo-protective agents, Hepes, and long-term preservatives (0.02% Sodium Azide). |
| Isotype             | IgG                                                                                                                                                                  |
| Storage Instruction | Store at -20°C for long term storage. Avoid freeze-thaw cycles.                                                                                                      |

### TARGET INFORMATION

|                  |                             |
|------------------|-----------------------------|
| Gene ID          | <a href="#">170643</a>      |
| Gene Symbol      | <a href="#">Kirrel1</a>     |
| UniProt ID       | <a href="#">KIRR1_MOUSE</a> |
| Immunogen Region | 200-250                     |
| Specificity      | Extracellular domain.       |

### ADDITIONAL INFORMATION

|      |                                                                                                                         |
|------|-------------------------------------------------------------------------------------------------------------------------|
| Note | STRICTLY FOR FURTHER SCIENTIFIC RESEARCH USE ONLY (RUO). MUST NOT TO BE USED IN DIAGNOSTIC OR THERAPEUTIC APPLICATIONS. |
|------|-------------------------------------------------------------------------------------------------------------------------|



Immunohistochemistry of mouse kidney with Anti-KIRREL primary antibody (STJ501885) at 1:100 dilution in Buffer. Section was treated with DAB and Haematoxylin stain. Magnification is at 40X.