

Anti-D4S234E antibody (140-190) {FITC} (STJ501974)

GENERAL INFORMATION

| | |
|---------------|--------------------|
| Product Type | Primary antibodies |
| Applications | ELISA/IHC/IP/WB |
| Host / Source | Rabbit |
| Reactivity | Human/Mouse/Rat |

PRODUCT PROPERTIES

| | |
|---------------------|--|
| Clonality | Polyclonal |
| Concentration | 0.5-1.5 $\mu\text{g}/\mu\text{l}$ |
| Conjugation | FITC |
| Purification | Affinity Purified |
| Dilution Range | WB: 1:500 ELISA: 1:10,000 IP: 1:200 IHC: 1:50-1:150 ICC: 1:50-1:150 |
| Formulation | Contains Tris, HCl/Glycine buffer pH 7.4-7.8, 30% Glycerol and 0.5% BSA, along with cryo-protective agents, Hepes, and long-term preservatives (0.02% Sodium Azide). |
| Isotype | IgG |
| Storage Instruction | Store at -20°C for long term storage. Avoid freeze-thaw cycles. |

TARGET INFORMATION

| | |
|------------------|--------------------------|
| Gene ID | 25247 |
| Gene Symbol | Nsg1 |
| UniProt ID | NSG1_RAT |
| Immunogen Region | 140-190 |

ADDITIONAL INFORMATION

| | |
|------|---|
| Note | STRICTLY FOR FURTHER SCIENTIFIC RESEARCH USE ONLY (RUO). MUST NOT TO BE USED IN DIAGNOSTIC OR THERAPEUTIC APPLICATIONS. |
|------|---|