

Anti-P2Y12 antibody (300-342) {Biotin} (STJ502042)

GENERAL INFORMATION		
	Product Type	Primary antibodies
	Applications	CM/ELISA/ICC/IF/IHC/IP/WB
	Host / Source	Rabbit
	Reactivity	Chimpanzee/Human/Monkey/Rat
PRODUCT PROPERTIES		
	Clonality	Polyclonal
	Concentration	0.57-0.68 µg/µl
	Conjugation	Biotin
	Purification	Affinity Purified
	Dilution Range	WB: 1:500 ELISA: 1:10, 000 IP: 1:200 IHC: 1:100 ICC: 1:100 IF: 1:100
	Formulation	Contains Tris, HCl/Glycine buffer pH 7.4-7.8, 30% Glycerol and 0.5% BSA, along with cryo-protective agents, Hepes, and long-term preservatives (0.02% Sodium Azide).
	Isotype	IgG
	Storage Instruction	Store at-20°C for long term storage. Avoid freeze-thaw cycles.
TARGET INFORMATION		
	Immunogen Region	300-342
ADDITIONAL INFORMATION		
	Note	STRICTLY FOR FURTHER SCIENTIFIC RESEARCH USE ONLY (RUO). MUST NOT TO BE USED IN DIAGNOSTIC OR THERAPEUTIC APPLICATIONS.