

**Anti-P2Y2 antibody (360-410) {Biotin} (STJ502045)****GENERAL INFORMATION**

|               |                           |
|---------------|---------------------------|
| Product Type  | Primary antibodies        |
| Applications  | CM/ELISA/ICC/IF/IHC/IP/WB |
| Host / Source | Rabbit                    |
| Reactivity    | Cat/Human/Mouse/Rat       |

**PRODUCT PROPERTIES**

|                     |  |
|---------------------|--|
| Clonality           | Polyclonal   |
| Concentration       | 0.52-0.68 $\mu\text{g}/\mu\text{l}$  |
| Conjugation         | Biotin   |
| Purification        | Affinity Purified  |
| Dilution Range      | WB: 1:500<br>ELISA: 1:10,000<br>IP: 1:200<br>IHC: 1:100<br>ICC: 1:100<br>IF: 1:100   |
| Formulation         | Contains Tris, HCl/Glycine buffer pH 7.4-7.8, 30% Glycerol and 0.5% BSA, along with cryo-protective agents, Hepes, and long-term preservatives (0.02% Sodium Azide). |
| Isotype             | IgG  |
| Storage Instruction | Store at -20°C for long term storage. Avoid freeze-thaw cycles.  |

**TARGET INFORMATION**

|                  |         |
|------------------|---------|
| Immunogen Region | 360-410 |
|------------------|---------|

**ADDITIONAL INFORMATION**

|      |   |
|------|---|
| Note | STRICTLY FOR FURTHER SCIENTIFIC RESEARCH USE ONLY (RUO). MUST NOT TO BE USED IN DIAGNOSTIC OR THERAPEUTIC APPLICATIONS. |
|------|---|