

Anti-Phospho-BTK antibody (544-561) {FITC} (STJ502178)

GENERAL INFORMATION

Product Type	Primary antibodies
Applications	ELISA/WB
Host / Source	Rabbit
Reactivity	Human/Monkey/Mouse/Pig/Rat

PRODUCT PROPERTIES

Clonality	Polyclonal
Concentration	0.6 µg/µl
Conjugation	FITC
Purification	Affinity Purified
Dilution Range	WB: 1:500 DB: 1:10, 000 ELISA: 1:10, 000
Formulation	Contains Tris, HCl/Glycine buffer pH 7.4-7.8, 30% Glycerol and 0.5% BSA, along with cryo-protective agents, Heps, and long-term preservatives (0.02% Sodium Azide).
Isotype	IgG
Storage Instruction	Store at -20°C for long term storage. Avoid freeze-thaw cycles.

TARGET INFORMATION

Gene ID	695
Gene Symbol	BTK
UniProt ID	BTK_HUMAN
Immunogen Region	544-561
Specificity	Does not cross react to any other cellular protein.

ADDITIONAL INFORMATION

Note	STRICTLY FOR FURTHER SCIENTIFIC RESEARCH USE ONLY (RUO). MUST NOT TO BE USED IN DIAGNOSTIC OR THERAPEUTIC APPLICATIONS.
------	---