

Anti-PDE1A antibody (51-101) {Biotin} (STJ502230)

GENERAL INFORMATION		
	Product Type	Primary antibodies
	Applications	CM/ELISA/ICC/IF/IHC/IP/WB
	Host / Source	Rabbit
	Reactivity	Human/Mouse/Rat
PRODUCT PROPERTIES		
	Clonality	Polyclonal
	Concentration	0.59 µg/µl
	Conjugation	Biotin
	Purification	Affinity Purified
	Dilution Range	WB: 1:500-1:1, 000 DB: 1:10, 000 ELISA: 1:10, 000 IP: 1:200 IHC: 1:200 ICC: 1:200 IF: 1:200 CFM: 1:200
	Formulation	Contains Tris, HCl/Glycine buffer pH 7.4-7.8, 30% Glycerol and 0.5% BSA, along with cryo-protective agents, Hepes, and long-term preservatives (0.02% Sodium Azide).
	Isotype	IgG
	Storage Instruction	Store at -20°C for long term storage. Avoid freeze-thaw cycles.
TARGET INFORMATION		
	Immunogen Region	51-101
	Specificity	This antibody will not label other PDE family members. This antibody will label all PDE1 variants including PDE1A1, PDE1A3, PDE1A5.
ADDITIONAL INFORMATION		
	Note	STRICTLY FOR FURTHER SCIENTIFIC RESEARCH USE ONLY (RUO). MUST NOT TO BE USED IN DIAGNOSTIC OR THERAPEUTIC APPLICATIONS.