

## Anti-PDE2A antibody (C-Term) {Biotin} (STJ502246)

### GENERAL INFORMATION

Product Type	Primary antibodies
Applications	CM/ELISA/FRET/ICC/IF/IHC/IP/WB
Host / Source	Rabbit
Reactivity	Human/Mouse/Rat

### PRODUCT PROPERTIES

Clonality	Polyclonal
Concentration	0.7 µg/µl
Conjugation	Biotin
Purification	Affinity Purified
Dilution Range	WB: 1:250-1:1, 000 DB: 1:10, 000 ELISA: 1:10, 000 IP: 1:200 IHC: 1:250-1:1, 000 ICC: 1:250-1:1, 000 IF: 1:150-1:1, 000 CFM: 1:250-1:1, 000
Formulation	Contains Tris, HCl/Glycine buffer pH 7.4-7.8, 30% Glycerol and 0.5% BSA, along with cryo-protective agents, HEPES, and long-term preservatives (0.02% Sodium Azide).
Isotype	IgG
Storage Instruction	Store at -20°C for long term storage. Avoid freeze-thaw cycles.

### TARGET INFORMATION

Immunogen Region	C-Term
Specificity	This antibody will label PDE2A1, PDE2A2, PDE2A3 and PDE2A4 variants.

### ADDITIONAL INFORMATION

Note	STRICTLY FOR FURTHER SCIENTIFIC RESEARCH USE ONLY (RUO). MUST NOT TO BE USED IN DIAGNOSTIC OR THERAPEUTIC APPLICATIONS.
------	---