

Anti-PDE2A antibody (C-Term) {FITC} (STJ502247)

GENERAL INFORMATION	
Product Type	Primary antibodies
Applications	CM/ELISA/FRET/ICC/IF/IHC/IP/WB
Host / Source	Rabbit
Reactivity	Human/Mouse/Rat
PRODUCT PROPERTIES	
Clonality	Polyclonal
Concentration	0.7 µg/µl
Conjugation	FITC
Purification	Affinity Purified
Dilution Range	WB: 1:250-1:1, 000 DB: 1:10, 000 ELISA: 1:10, 000 IP: 1:200 IHC: 1:250-1:1, 000 ICC: 1:250-1:1, 000 IF: 1:150-1:1, 000 CFM: 1:250-1:1, 000
Formulation	Contains Tris, HCl/Glycine buffer pH 7.4-7.8, 30% Glycerol and 0.5% BSA, along with cryo-protective agents, Hepes, and long-term preservatives (0.02% Sodium Azide).
Isotype	IgG
Storage Instruction	Store at -20°C for long term storage. Avoid freeze-thaw cycles.
TARGET INFORMATION	
Immunogen Region	C-Term
Specificity	This antibody will label PDE2A1, PDE2A2, PDE2A3 and PDE2A4 variants.
ADDITIONAL INFORMATION	
Note	STRICTLY FOR FURTHER SCIENTIFIC RESEARCH USE ONLY (RUO). MUST NOT TO BE USED IN DIAGNOSTIC OR THERAPEUTIC APPLICATIONS.