

Anti-PDE3B antibody (1060-1110) (STJ502269)

GENERAL INFORMATION

Product Type	Primary antibodies
Applications	ELISA/IP/WB
Host / Source	Goat
Reactivity	Chicken/Human/Monkey/Mouse/Rat

PRODUCT PROPERTIES


Clonality	Polyclonal
Concentration	0.55 µg/µl
Conjugation	Unconjugated
Purification	Affinity Purified
Dilution Range	WB: 1:500 DB: 1:10, 000 ELISA: 1:10, 000 IP: 1:200 IHC: 1:500 ICC: 1:500 IF: 1:500 CFM: 1:500
Formulation	Contains Tris, HCl/Glycine buffer pH 7.4-7.8, 30% Glycerol and 0.5% BSA, along with cryo-protective agents, Hepes, and long-term preservatives (0.02% Sodium Azide).
Isotype	IgG
Storage Instruction	Store at -20°C for long term storage. Avoid freeze-thaw cycles.

TARGET INFORMATION

Gene ID	29516
Gene Symbol	Pde3b
UniProt ID	PDE3B_RAT
Immunogen Region	1060-1110
Specificity	Does not cross react with PDE4A variants nor with any other PDE family members. This antibody will recognize PDE3B variants.

ADDITIONAL INFORMATION

Note	STRICTLY FOR FURTHER SCIENTIFIC RESEARCH USE ONLY (RUO). MUST NOT TO BE USED IN DIAGNOSTIC OR THERAPEUTIC APPLICATIONS.
------	---



Western blot with Anti-PDE3B Antibody (STJ502269) on PD38 peptide. Antibody at 1:500 dilution in buffer.