

Anti-PDE5A antibody {FITC} (STJ502355)

GENERAL INFORMATION	
Product Type	Primary antibodies
Applications	CM/ELISA/ICC/IF/IHC/IP/WB
Host / Source	Rabbit
Reactivity	Chicken/Human
PRODUCT PROPERTIES	
Clonality	Polyclonal
Concentration	0.55-0.60 µg/µl
Conjugation	FITC
Purification	Affinity Purified
Dilution Range	WB: 1:500-1:1, 000 DB: 1:10, 000 ELISA: 1:10, 000 IP: 1:250 IHC: 1:250 ICC: 1:250 IF: 1:250 CFM: 1:250
Formulation	Contains Tris, HCl/Glycine buffer pH 7.4-7.8, 30% Glycerol and 0.5% BSA, along with cryo-protective agents, Hepes, and long-term preservatives (0.02% Sodium Azide).
Isotype	IgG
Storage Instruction	Store at -20°C for long term storage. Avoid freeze-thaw cycles.
TARGET INFORMATION	
Gene ID	8654
Gene Symbol	PDE5A
UniProt ID	PDE5A_HUMAN
ADDITIONAL INFORMATION	
Note	STRICTLY FOR FURTHER SCIENTIFIC RESEARCH USE ONLY (RUO). MUST NOT TO BE USED IN DIAGNOSTIC OR THERAPEUTIC APPLICATIONS.