

Anti-PDE9A antibody {Biotin} (STJ502390)

GENERAL INFORMATION		
	Product Type	Primary antibodies
	Applications	CM/ELISA/ICC/IF/IHC/IP/WB
	Host / Source	Rabbit
	Reactivity	Human/Mouse/Rat
PRODUCT PROPERTIES		
	Clonality	Polyclonal
	Concentration	0.55 µg/µl
	Conjugation	Biotin
	Purification	Affinity Purified
	Dilution Range	WB: 1:500 DB: 1:100 ELISA: 1:100 IP: 1:250 IHC: 1:100 ICC: 1:100 IF: 1:100 CFM: 1:100
	Formulation	Contains Tris, HCl/Glycine buffer pH 7.4-7.8, 30% Glycerol and 0.5% BSA, along with cryo-protective agents, Hepes, and long-term preservatives (0.02% Sodium Azide).
	Isotype	IgG
	Storage Instruction	Store at-20°C for long term storage. Avoid freeze-thaw cycles.
TARGET INFORMATION		
	Gene ID	5152
	Gene Symbol	PDE9A
	UniProt ID	PDE9A_HUMAN
	Specificity	Does not cross-react with PDE9B protein or any other PDE family members.
ADDITIONAL INFORMATION		
	Note	STRICTLY FOR FURTHER SCIENTIFIC RESEARCH USE ONLY (RUO). MUST NOT TO BE USED IN DIAGNOSTIC OR THERAPEUTIC APPLICATIONS.