

## Anti-PDE6C antibody (1-75) {Biotin} (STJ502396)

| GENERAL INFORMATION    |                                                                                                                                                                      |
|------------------------|----------------------------------------------------------------------------------------------------------------------------------------------------------------------|
| Product Type           | Primary antibodies                                                                                                                                                   |
| Applications           | ELISA/IP/WB                                                                                                                                                          |
| Host / Source          | Rabbit                                                                                                                                                               |
| Reactivity             | Bovine/Human/Mouse/Rat                                                                                                                                               |
| PRODUCT PROPERTIES     |                                                                                                                                                                      |
| Clonality              | Polyclonal                                                                                                                                                           |
| Concentration          | 0.50-0.75 µg/µl                                                                                                                                                      |
| Conjugation            | Biotin                                                                                                                                                               |
| Purification           | Affinity Purified                                                                                                                                                    |
| Dilution Range         | WB: 1:500<br>DB: 1:10, 000<br>ELISA: 1:10, 000<br>IP: 1:200<br>IHC: 1:50-1:200                                                                                       |
| Formulation            | Contains Tris, HCl/Glycine buffer pH 7.4-7.8, 30% Glycerol and 0.5% BSA, along with cryo-protective agents, Hepes, and long-term preservatives (0.02% Sodium Azide). |
| Isotype                | IgG                                                                                                                                                                  |
| Storage Instruction    | Store at -20°C for long term storage. Avoid freeze-thaw cycles.                                                                                                      |
| TARGET INFORMATION     |                                                                                                                                                                      |
| Gene ID                | <a href="#">5146</a>                                                                                                                                                 |
| Gene Symbol            | <a href="#">PDE6C</a>                                                                                                                                                |
| UniProt ID             | <a href="#">PDE6C_HUMAN</a>                                                                                                                                          |
| Immunogen Region       | 1-75                                                                                                                                                                 |
| ADDITIONAL INFORMATION |                                                                                                                                                                      |
| Note                   | STRICTLY FOR FURTHER SCIENTIFIC RESEARCH USE ONLY (RUO). MUST NOT TO BE USED IN DIAGNOSTIC OR THERAPEUTIC APPLICATIONS.                                              |