

Anti-PMS2 antibody (220-270) {Biotin} (STJ502563)

GENERAL INFORMATION		
	Product Type	Primary antibodies
	Applications	ELISA/WB
	Host / Source	Rabbit
	Reactivity	Human/Mouse/Rat
PRODUCT PROPERTIES		
	Clonality	Polyclonal
	Concentration	0.65-0.69 µg/µl
	Conjugation	Biotin
	Purification	Affinity Purified
	Dilution Range	WB: 1:500-1:1, 000 DB: 1:20, 000-1:50, 000 ELISA: 1:20, 000-1:50, 000 IP: 1:50-1:200 IHC: 1:50-1:200
	Formulation	Contains Tris, HCl/Glycine buffer pH 7.4-7.8, 30% Glycerol and 0.5% BSA, along with cryo-protective agents, Hepes, and long-term preservatives (0.02% Sodium Azide).
	Isotype	IgG
	Storage Instruction	Store at-20°C for long term storage. Avoid freeze-thaw cycles.
TARGET INFORMATION		
	Gene ID	5395
	Gene Symbol	PMS2
	UniProt ID	PMS2_HUMAN
	Immunogen Region	220-270
ADDITIONAL INFORMATION		
	Note	STRICTLY FOR FURTHER SCIENTIFIC RESEARCH USE ONLY (RUO). MUST NOT TO BE USED IN DIAGNOSTIC OR THERAPEUTIC APPLICATIONS.