

Anti-Phospho-PDE3B antibody (260-310) {Biotin} (STJ502610)

GENERAL INFORMATION		
	Product Type	Primary antibodies
	Applications	ELISA/IP/WB
	Host / Source	Rabbit
	Reactivity	Chicken/Human/Monkey/Mouse/Rat
PRODUCT PROPERTIES		
	Clonality	Polyclonal
	Concentration	0.70 µg/µl
	Conjugation	Biotin
	Purification	Affinity Purified
	Dilution Range	WB: 1:500 DB: 1:10, 000 ELISA: 1:10, 000
	Formulation	Contains Tris, HCl/Glycine buffer pH 7.4-7.8, 30% Glycerol and 0.5% BSA, along with cryo-protective agents, Hepes, and long-term preservatives (0.02% Sodium Azide).
	Isotype	IgG
	Storage Instruction	Store at-20°C for long term storage. Avoid freeze-thaw cycles.
TARGET INFORMATION		
	Gene ID	5140
	Gene Symbol	PDE3B
	UniProt ID	PDE3B_HUMAN
	Immunogen Region	260-310
ADDITIONAL INFORMATION		
	Note	STRICTLY FOR FURTHER SCIENTIFIC RESEARCH USE ONLY (RUO). MUST NOT TO BE USED IN DIAGNOSTIC OR THERAPEUTIC APPLICATIONS.