

Anti-Phospho-PDE5A antibody (120-170) {Biotin} (STJ502622)

GENERAL INFORMATION		
	Product Type	Primary antibodies
	Applications	ELISA/ICC/IHC/IP/WB
	Host / Source	Rabbit
	Reactivity	Bovine/Human/Mouse/Pig/Rat/Sheep
PRODUCT PROPERTIES		
	Clonality	Polyclonal
	Concentration	0.50 µg/µl
	Conjugation	Biotin
	Purification	Affinity Purified
	Dilution Range	WB: 1:500 DB: 1:10, 000 ELISA: 1:10, 000 IP: 1:200 IHC: 1:100 ICC: 1:100 IF: 1:100 CFM: 1:100
	Formulation	Contains Tris, HCl/Glycine buffer pH 7.4-7.8, 30% Glycerol and 0.5% BSA, along with cryo-protective agents, Hepes, and long-term preservatives (0.02% Sodium Azide).
	Isotype	IgG
	Storage Instruction	Store at-20°C for long term storage. Avoid freeze-thaw cycles.
TARGET INFORMATION		
	Immunogen Region	120-170
	Specificity	This antibody will detect only PDE5A protein that has been phosphorylated at the PKG site. Will not detect unphosphorylated PDE5A, or other PDE family members.
ADDITIONAL INFORMATION		
	Note	STRICTLY FOR FURTHER SCIENTIFIC RESEARCH USE ONLY (RUO). MUST NOT TO BE USED IN DIAGNOSTIC OR THERAPEUTIC APPLICATIONS.