

## Anti-PTPD1 antibody (400-500) (STJ502717)

### GENERAL INFORMATION

|               |                                |
|---------------|--------------------------------|
| Product Type  | Primary antibodies             |
| Applications  | ELISA/IHC/IP/WB                |
| Host / Source | Rabbit                         |
| Reactivity    | Chicken/Human/Monkey/Mouse/Rat |

### PRODUCT PROPERTIES


|                     |                                                                                                                                                                      |
|---------------------|----------------------------------------------------------------------------------------------------------------------------------------------------------------------|
| Clonality           | Polyclonal                                                                                                                                                           |
| Concentration       | 0.65-0.87 $\mu\text{g}/\mu\text{l}$                                                                                                                                  |
| Conjugation         | Unconjugated                                                                                                                                                         |
| Purification        | Affinity Purified                                                                                                                                                    |
| Dilution Range      | WB: 1:250<br>DB: 1:10, 000<br>ELISA: 1:10, 000<br>IP: 1:200<br>IHC: 1:50-1:150                                                                                       |
| Formulation         | Contains Tris, HCl/Glycine buffer pH 7.4-7.8, 30% Glycerol and 0.5% BSA, along with cryo-protective agents, Hepes, and long-term preservatives (0.02% Sodium Azide). |
| Isotype             | IgG                                                                                                                                                                  |
| Storage Instruction | Store at -20°C for long term storage. Avoid freeze-thaw cycles.                                                                                                      |

### TARGET INFORMATION

|                  |         |
|------------------|---------|
| Immunogen Region | 400-500 |
|------------------|---------|

### ADDITIONAL INFORMATION

|      |                                                                                                                         |
|------|-------------------------------------------------------------------------------------------------------------------------|
| Note | STRICTLY FOR FURTHER SCIENTIFIC RESEARCH USE ONLY (RUO). MUST NOT TO BE USED IN DIAGNOSTIC OR THERAPEUTIC APPLICATIONS. |
|------|-------------------------------------------------------------------------------------------------------------------------|



Immunohistochemistry of mouse spleen with Anti-PTPD1 primary antibody (STJ502717) at 1:50 dilution in Buffer. Section was treated with DAB and Haematoxylin stain. Magnification is at 40X.