

## Anti-ROCK1 antibody (850-900) {Biotin} (STJ502820)

### GENERAL INFORMATION

|               |                            |
|---------------|----------------------------|
| Product Type  | Primary antibodies         |
| Applications  | ELISA/IHC/IP/WB            |
| Host / Source | Rabbit                     |
| Reactivity    | Human/Monkey/Mouse/Pig/Rat |

### PRODUCT PROPERTIES

|                     |  |
|---------------------|--|
| Clonality           | Polyclonal   |
| Concentration       | 0.57 µg/µl   |
| Conjugation         | Biotin   |
| Purification        | Affinity Purified  |
| Dilution Range      | WB: 1:500<br>DB: 1:10, 000<br>ELISA: 1:10, 000<br>IP: 1:200<br>IHC: 1:50-1:100   |
| Formulation         | Contains Tris, HCl/Glycine buffer pH 7.4-7.8, 30% Glycerol and 0.5% BSA, along with cryo-protective agents, Hepes, and long-term preservatives (0.02% Sodium Azide). |
| Isotype             | IgG  |
| Storage Instruction | Store at -20°C for long term storage. Avoid freeze-thaw cycles.  |

### TARGET INFORMATION

|                  |                             |
|------------------|-----------------------------|
| Gene ID          | <a href="#">6093</a>        |
| Gene Symbol      | <a href="#">ROCK1</a>       |
| UniProt ID       | <a href="#">ROCK1_HUMAN</a> |
| Immunogen Region | 850-900                     |

### ADDITIONAL INFORMATION

|      |   |
|------|---|
| Note | STRICTLY FOR FURTHER SCIENTIFIC RESEARCH USE ONLY (RUO). MUST NOT TO BE USED IN DIAGNOSTIC OR THERAPEUTIC APPLICATIONS. |
|------|---|