

**Anti-C8orf44-SGK3 antibody (440-490) {FITC} (STJ502986)****GENERAL INFORMATION**

Product Type	Primary antibodies
Applications	CM/ELISA/IC/IF/IHC/IP/WB
Host / Source	Chicken
Reactivity	Human/Mouse/Rat

**PRODUCT PROPERTIES**

Clonality	Polyclonal
Concentration	0.6 µg/µl
Conjugation	FITC
Purification	Affinity Purified
Dilution Range	WB: 1:500 DB: 1:10, 000 ELISA: 1:10, 000 IP: 1:200 IHC: 1:100 ICC: 1:100 IF: 1:100 CFM: 1:100
Formulation	Contains Tris, HCl/Glycine buffer pH 7.4-7.8, 30% Glycerol and 0.5% BSA, along with cryo-protective agents, HEPES, and long-term preservatives (0.02% Sodium Azide).
Isotype	IgY
Storage Instruction	Store at -20°C for long term storage. Avoid freeze-thaw cycles.

**TARGET INFORMATION**

Gene ID	<a href="#">100533105/23678</a>
Gene Symbol	<a href="#">SGK3</a>
UniProt ID	<a href="#">SGK3_HUMAN</a>
Immunogen Region	440-490

**ADDITIONAL INFORMATION**

Note	STRICTLY FOR FURTHER SCIENTIFIC RESEARCH USE ONLY (RUO). MUST NOT TO BE USED IN DIAGNOSTIC OR THERAPEUTIC APPLICATIONS.
------	---