

Anti-Syncytin 2 antibody (380-430) (STJ503173)

GENERAL INFORMATION

Product Type	Primary antibodies
Applications	ELISA/IP/WB
Host / Source	Rabbit
Reactivity	Mouse/Rat

PRODUCT PROPERTIES


Clonality	Polyclonal
Concentration	0.5 µg/µl
Conjugation	Unconjugated
Purification	Affinity Purified
Dilution Range	WB: 1:500 ELISA: 1:50, 000 IP: 1:200
Formulation	Contains Tris, HCl/Glycine buffer pH 7.4-7.8, 30% Glycerol and 0.5% BSA, along with cryo-protective agents, Hepes, and long-term preservatives (0.02% Sodium Azide).
Isotype	IgG
Storage Instruction	Store at -20°C for long term storage. Avoid freeze-thaw cycles.

TARGET INFORMATION

Gene ID	239167
Gene Symbol	Synb
UniProt ID	SYNB_MOUSE
Immunogen Region	380-430

ADDITIONAL INFORMATION

Note STRICTLY FOR FURTHER SCIENTIFIC RESEARCH USE ONLY (RUO). MUST NOT TO BE USED IN DIAGNOSTIC OR THERAPEUTIC APPLICATIONS.



Western blot with Anti-Syncytin 2 Antibody (STJ503173) on rat: 1) Cortex, 2) Midbrain, 3) Cerebellum, 4) Hippocampus, 5) Kidney, 6) Adrenal gland, 7) Thymus, 8) Trachea, 9) Ovary, 10) Bone marrow, 11) Testis. Antibody at 1:250 dilution in buffer.