

Anti-SYDE1 antibody (680-735) {Biotin} (STJ503180)

GENERAL INFORMATION	
Product Type	Primary antibodies
Applications	ELISA/WB
Host / Source	Rabbit
Reactivity	C.elegans/Human/Mouse
PRODUCT PROPERTIES	
Clonality	Polyclonal
Concentration	0.67-0.68 µg/µl
Conjugation	Biotin
Purification	Affinity Purified
Dilution Range	WB: 1:500 DB: 1:10, 000 ELISA: 1:10, 000 IP: 1:200 IHC: 1:50-1:150
Formulation	Contains Tris, HCl/Glycine buffer pH 7.4-7.8, 30% Glycerol and 0.5% BSA, along with cryo-protective agents, Hepes, and long-term preservatives (0.02% Sodium Azide).
Isotype	IgG
Storage Instruction	Store at -20°C for long term storage. Avoid freeze-thaw cycles.
TARGET INFORMATION	
Gene ID	85360
Gene Symbol	SYDE1
UniProt ID	SYDE1_HUMAN
Immunogen Region	680-735
ADDITIONAL INFORMATION	
Note	STRICTLY FOR FURTHER SCIENTIFIC RESEARCH USE ONLY (RUO). MUST NOT TO BE USED IN DIAGNOSTIC OR THERAPEUTIC APPLICATIONS.