

Anti-TGM1 antibody {FITC} (STJ503242)

GENERAL INFORMATION		
	Product Type	Primary antibodies
	Applications	CM/ELISA/ICC/IHC/IF/WB
	Host / Source	Rabbit
	Reactivity	Human/Mouse/Rat
PRODUCT PROPERTIES		
	Clonality	Polyclonal
	Concentration	0.61-0.65 µg/µl
	Conjugation	FITC
	Purification	Affinity Purified
	Dilution Range	WB: 1:450 DB: 1:30, 000 ELISA: 1:30, 000 IHC: 1:250 ICC: 1:250 IF: 1:250 CFM: 1:250
	Formulation	Contains Tris, HCl/Glycine buffer pH 7.4-7.8, 30% Glycerol and 0.5% BSA, along with cryo-protective agents, Hepes, and long-term preservatives (0.02% Sodium Azide).
	Isotype	IgG
	Storage Instruction	Store at-20°C for long term storage. Avoid freeze-thaw cycles.
TARGET INFORMATION		
	Gene ID	74145
	Gene Symbol	F13a1
	UniProt ID	F13A_MOUSE
ADDITIONAL INFORMATION		
	Note	STRICTLY FOR FURTHER SCIENTIFIC RESEARCH USE ONLY (RUO). MUST NOT TO BE USED IN DIAGNOSTIC OR THERAPEUTIC APPLICATIONS.