

Anti-TEK antibody (540-590) {Biotin} (STJ503261)

GENERAL INFORMATION

| | |
|---------------|--------------------|
| Product Type | Primary antibodies |
| Applications | DB/ELISA/IHC/IP/WB |
| Host / Source | Chicken |
| Reactivity | Mouse/Rat |

PRODUCT PROPERTIES

| | |
|---------------------|--|
| Clonality | Polyclonal |
| Concentration | 0.68-0.70 $\mu\text{g}/\mu\text{l}$ |
| Conjugation | Biotin |
| Purification | Affinity Purified |
| Dilution Range | WB: 1:500 DB: 1:4, 000 ELISA: 1:4, 000 IP: 1:150 IHC: 1:50-1:100 |
| Formulation | Contains Tris, HCl/Glycine buffer pH 7.4-7.8, 30% Glycerol and 0.5% BSA, along with cryo-protective agents, Hepes, and long-term preservatives (0.02% Sodium Azide). |
| Isotype | IgY |
| Storage Instruction | Store at -20°C for long term storage. Avoid freeze-thaw cycles. |

TARGET INFORMATION

| | |
|------------------|----------------------------|
| Gene ID | 7010 |
| Gene Symbol | TEK |
| UniProt ID | TIE2_HUMAN |
| Immunogen Region | 540-590 |

ADDITIONAL INFORMATION

| | |
|------|---|
| Note | STRICTLY FOR FURTHER SCIENTIFIC RESEARCH USE ONLY (RUO). MUST NOT TO BE USED IN DIAGNOSTIC OR THERAPEUTIC APPLICATIONS. |
|------|---|