

Anti-TICAM2 antibody (140-170) {Biotin} (STJ503267)

GENERAL INFORMATION

Product Type	Primary antibodies
Applications	ELISA/WB
Host / Source	Rabbit
Reactivity	Cat/Human/Monkey/Pig

PRODUCT PROPERTIES

Clonality	Polyclonal
Concentration	0.60-0.65 $\mu\text{g}/\mu\text{l}$
Conjugation	Biotin
Purification	Affinity Purified
Dilution Range	WB: 1:250 DB: 1:10, 000 ELISA: 1:10, 000 IP: 1:200
Formulation	Contains Tris, HCl/Glycine buffer pH 7.4-7.8, 30% Glycerol and 0.5% BSA, along with cryo-protective agents, Hepes, and long-term preservatives (0.02% Sodium Azide).
Isotype	IgG
Storage Instruction	Store at -20°C for long term storage. Avoid freeze-thaw cycles.

TARGET INFORMATION

Gene ID	100302736/353376
Gene Symbol	TICAM2
UniProt ID	TCAM2_HUMAN
Immunogen Region	140-170

ADDITIONAL INFORMATION

Note	STRICTLY FOR FURTHER SCIENTIFIC RESEARCH USE ONLY (RUO). MUST NOT TO BE USED IN DIAGNOSTIC OR THERAPEUTIC APPLICATIONS.
------	---