

## Anti-TOMM6 antibody (24-74) {FITC} (STJ503353)

### GENERAL INFORMATION

|               |                               |
|---------------|-------------------------------|
| Product Type  | Primary antibodies            |
| Applications  | ELISA/IP/WB                   |
| Host / Source | Rabbit                        |
| Reactivity    | Bovine/Human/Monkey/Mouse/Rat |

### PRODUCT PROPERTIES

|                     |                                                                                                                                                                      |
|---------------------|----------------------------------------------------------------------------------------------------------------------------------------------------------------------|
| Clonality           | Polyclonal                                                                                                                                                           |
| Concentration       | 0.55-0.75 µg/µl                                                                                                                                                      |
| Conjugation         | FITC                                                                                                                                                                 |
| Purification        | Affinity Purified                                                                                                                                                    |
| Dilution Range      | WB: 1:500<br>DB: 1:10, 000<br>ELISA: 1:10, 000<br>IP: 1:150                                                                                                          |
| Formulation         | Contains Tris, HCl/Glycine buffer pH 7.4-7.8, 30% Glycerol and 0.5% BSA, along with cryo-protective agents, Hepes, and long-term preservatives (0.02% Sodium Azide). |
| Isotype             | IgG                                                                                                                                                                  |
| Storage Instruction | Store at -20°C for long term storage. Avoid freeze-thaw cycles.                                                                                                      |

### TARGET INFORMATION

|                  |                            |
|------------------|----------------------------|
| Gene ID          | <a href="#">100188893</a>  |
| Gene Symbol      | <a href="#">TOMM6</a>      |
| UniProt ID       | <a href="#">TOM6_HUMAN</a> |
| Immunogen Region | 24-74                      |

### ADDITIONAL INFORMATION

|      |                                                                                                                         |
|------|-------------------------------------------------------------------------------------------------------------------------|
| Note | STRICTLY FOR FURTHER SCIENTIFIC RESEARCH USE ONLY (RUO). MUST NOT TO BE USED IN DIAGNOSTIC OR THERAPEUTIC APPLICATIONS. |
|------|-------------------------------------------------------------------------------------------------------------------------|