

Anti-ULK1 antibody (1000-1050) {Biotin} (STJ503462)

GENERAL INFORMATION

Product Type	Primary antibodies
Applications	ELISA/WB
Host / Source	Rabbit
Reactivity	Bovine/Human/Monkey/Mouse/Panda

PRODUCT PROPERTIES

Clonality	Polyclonal
Concentration	0.58 µg/µl
Conjugation	Biotin
Purification	Affinity Purified
Dilution Range	WB: 1:250 ELISA: 1:10, 000 IP: 1:200 IF: 1:100
Formulation	Contains Tris, HCl/Glycine buffer pH 7.4-7.8, 30% Glycerol and 0.5% BSA, along with cryo-protective agents, Hepes, and long-term preservatives (0.02% Sodium Azide).
Isotype	IgG
Storage Instruction	Store at -20°C for long term storage. Avoid freeze-thaw cycles.

TARGET INFORMATION

Gene ID	8408
Gene Symbol	ULK1
UniProt ID	ULK1_HUMAN
Immunogen Region	1000-1050

ADDITIONAL INFORMATION

Note	STRICTLY FOR FURTHER SCIENTIFIC RESEARCH USE ONLY (RUO). MUST NOT TO BE USED IN DIAGNOSTIC OR THERAPEUTIC APPLICATIONS.
------	---