

## 5HT1B Blocking Peptide for STJ50004 peptide (STJ503613)

### GENERAL INFORMATION

|              |                                   |
|--------------|-----------------------------------|
| Product Type | Biomolecules                      |
| Applications | Immunodepletion/Immunocompetition |

### PRODUCT PROPERTIES

|                     |  |
|---------------------|--|
| Formulation         | Liquid form at 2.5mg/ml concentration in PBS. Up to 5% DMSO can be added. Orders with >1mg can be supplied in lyophilized powder form, or in buffer of choice. |
| Storage Instruction | Store at -20°C for long term storage. Avoid freeze-thaw cycles.  |

### TARGET INFORMATION

|                  |   |
|------------------|---|
| Gene ID          | <a href="#">15551</a>   |
| Gene Symbol      | <a href="#">Htr1b</a>   |
| UniProt ID       | <a href="#">5HT1B_MOUSE</a>   |
| Immunogen Region | 160-260   |
| Specificity      | This blocking peptide is recommended for use in combination with 5HT1B antibody, STJ50004 |

### ADDITIONAL INFORMATION

|      |   |
|------|---|
| Note | STRICTLY FOR FURTHER SCIENTIFIC RESEARCH USE ONLY (RUO). MUST NOT TO BE USED IN DIAGNOSTIC OR THERAPEUTIC APPLICATIONS. |
|------|---|