

## AKR1B1 Positive Control for STJ500068 peptide (STJ503640)

### GENERAL INFORMATION

Product Type	Biomolecules
Applications	WB

### PRODUCT PROPERTIES

Dilution Range	WB: 1:500
Formulation	Provided as 100 uL ready-to-use, in SDS-PAGE sample buffer (Laemelli's buffer) containing Tris, pH 6.8, 1 % SDS, Glycerol and Bromophenolblue blue as tracking dye. The sample is reduced by adding 2% beta mercaptoethanol. The protein concentration is
Storage Instruction	Store at-20°C for long term storage. Avoid freeze-thaw cycles.

### TARGET INFORMATION

Gene ID	<a href="#">231</a>
Gene Symbol	<a href="#">AKR1B1</a>
UniProt ID	<a href="#">ALDR_HUMAN</a>
Specificity	This is positive control is recommended for use in combination with AKR1B1 antibody STJ500068.

### ADDITIONAL INFORMATION

Note	STRICTLY FOR FURTHER SCIENTIFIC RESEARCH USE ONLY (RUO). MUST NOT TO BE USED IN DIAGNOSTIC OR THERAPEUTIC APPLICATIONS.
------	-------------------------------------------------------------------------------------------------------------------------