

ANGPTL5 Blocking Peptide for STJ500086 peptide (STJ503644)

GENERAL INFORMATION

Product Type	Biomolecules
Applications	Immunodepletion/Immunocompetition

PRODUCT PROPERTIES

Formulation	Liquid form at 2.5mg/ml concentration in PBS. Up to 5% DMSO can be added. Orders with >1mg can be supplied in lyophilized powder form, or in buffer of choice.
Storage Instruction	Store at -20°C for long term storage. Avoid freeze-thaw cycles.

TARGET INFORMATION

Gene ID	253935
Gene Symbol	ANGPTL5
UniProt ID	ANGL5_HUMAN
Immunogen Region	340-388
Specificity	This blocking peptide is recommended for use in combination with ANGPTL5 antibody, STJ500086

ADDITIONAL INFORMATION

Note	STRICTLY FOR FURTHER SCIENTIFIC RESEARCH USE ONLY (RUO). MUST NOT TO BE USED IN DIAGNOSTIC OR THERAPEUTIC APPLICATIONS.
------	---