

## BCKDH Blocking Peptide peptide (STJ503694)

### GENERAL INFORMATION

|              |                                   |
|--------------|-----------------------------------|
| Product Type | Biomolecules                      |
| Applications | Immunodepletion/Immunocompetition |

### PRODUCT PROPERTIES

|                     |                                                                                                                                                                |
|---------------------|----------------------------------------------------------------------------------------------------------------------------------------------------------------|
| Formulation         | Liquid form at 2.5mg/ml concentration in PBS. Up to 5% DMSO can be added. Orders with >1mg can be supplied in lyophilized powder form, or in buffer of choice. |
| Storage Instruction | Store at -20°C for long term storage. Avoid freeze-thaw cycles.                                                                                                |

### TARGET INFORMATION

|                  |                            |
|------------------|----------------------------|
| Gene ID          | <a href="#">10295</a>      |
| Gene Symbol      | <a href="#">BCKDK</a>      |
| UniProt ID       | <a href="#">BCKD_HUMAN</a> |
| Immunogen Region | 170-200                    |

### ADDITIONAL INFORMATION

|      |                                                                                                                         |
|------|-------------------------------------------------------------------------------------------------------------------------|
| Note | STRICTLY FOR FURTHER SCIENTIFIC RESEARCH USE ONLY (RUO). MUST NOT TO BE USED IN DIAGNOSTIC OR THERAPEUTIC APPLICATIONS. |
|------|-------------------------------------------------------------------------------------------------------------------------|