

## SLC7A10 Blocking Peptide for STJ500165 peptide (STJ503697)

### GENERAL INFORMATION

Product Type	Biomolecules
Applications	Immunodepletion/Immunocompetition

### PRODUCT PROPERTIES

Formulation	Liquid form at 2.5mg/ml concentration in PBS. Up to 5% DMSO can be added. Orders with >1mg can be supplied in lyophilized powder form, or in buffer of choice.
Storage Instruction	Store at -20°C for long term storage. Avoid freeze-thaw cycles.

### TARGET INFORMATION

Gene ID	<a href="#">53896</a>
Gene Symbol	<a href="#">Slc7a10</a>
UniProt ID	<a href="#">AAA1_MOUSE</a>
Immunogen Region	1-50
Specificity	This blocking peptide is recommended for use in combination with SLC7A10 antibody, STJ500165

### ADDITIONAL INFORMATION

Note	STRICTLY FOR FURTHER SCIENTIFIC RESEARCH USE ONLY (RUO). MUST NOT TO BE USED IN DIAGNOSTIC OR THERAPEUTIC APPLICATIONS.
------	---