

CD163 Blocking Peptide for STJ500423 peptide (STJ503789)

GENERAL INFORMATION	
Product Type	Biomolecules
Applications	ELISA/CM/ICC/IF/IHC/IP/WB
PRODUCT PROPERTIES	
Formulation	Liquid form at 2.5mg/ml concentration in PBS. Up to 5% DMSO can be added. Orders with >1mg can be supplied in lyophilized powder form, or in buffer of choice.
Storage Instruction	Store at -20°C for long term storage. Avoid freeze-thaw cycles.
TARGET INFORMATION	
Immunogen Region	600-650
Specificity	This blocking peptide is recommended for use in combination with CD163 antibody, STJ500423
ADDITIONAL INFORMATION	
Note	STRICTLY FOR FURTHER SCIENTIFIC RESEARCH USE ONLY (RUO). MUST NOT TO BE USED IN DIAGNOSTIC OR THERAPEUTIC APPLICATIONS.