

CaMKII gamma Positive Control for STJ500360 peptide (STJ503812)

GENERAL INFORMATION

| | |
|--------------|--------------|
| Product Type | Biomolecules |
| Applications | WB |

PRODUCT PROPERTIES

| | |
|---------------------|---|
| Dilution Range | WB: 1:500 |
| Formulation | Provided as 100 uL ready-to-use, in SDS-PAGE sample buffer (Laemelli's buffer) containing Tris, pH 6.8, 1 % SDS, Glycerol and Bromophenolblue blue as tracking dye. The sample is reduced by adding 2% beta mercaptoethanol. The protein concentration is |
| Storage Instruction | Store at-20°C for long term storage. Avoid freeze-thaw cycles. |

TARGET INFORMATION

| | |
|-------------|--|
| Specificity | This is positive control is recommended for use in combination with CaMKII gamma antibody STJ500360. |
|-------------|--|

ADDITIONAL INFORMATION

| | |
|------|---|
| Note | STRICTLY FOR FURTHER SCIENTIFIC RESEARCH USE ONLY (RUO). MUST NOT TO BE USED IN DIAGNOSTIC OR THERAPEUTIC APPLICATIONS. |
|------|---|