

NT5C2 Blocking Peptide for STJ500602 peptide (STJ503868)

GENERAL INFORMATION

Product Type	Biomolecules
Applications	Immunodepletion/Immunocompetition

PRODUCT PROPERTIES

Formulation	Liquid form at 2.5mg/ml concentration in PBS. Up to 5% DMSO can be added. Orders with >1mg can be supplied in lyophilized powder form, or in buffer of choice.
Storage Instruction	Store at -20°C for long term storage. Avoid freeze-thaw cycles.

TARGET INFORMATION

Gene ID	22978
Gene Symbol	NT5C2
UniProt ID	5NTC_HUMAN
Specificity	This blocking peptide is recommended for use in combination with NT5C2 antibody, STJ500602

ADDITIONAL INFORMATION

Note	STRICTLY FOR FURTHER SCIENTIFIC RESEARCH USE ONLY (RUO). MUST NOT TO BE USED IN DIAGNOSTIC OR THERAPEUTIC APPLICATIONS.
------	---