

## CXCR1 Blocking Peptide for STJ500647 peptide (STJ503880)

### GENERAL INFORMATION

|              |                                   |
|--------------|-----------------------------------|
| Product Type | Biomolecules                      |
| Applications | Immunodepletion/Immunocompetition |

### PRODUCT PROPERTIES

|                     |  |
|---------------------|--|
| Formulation         | Liquid form at 2.5mg/ml concentration in PBS. Up to 5% DMSO can be added. Orders with >1mg can be supplied in lyophilized powder form, or in buffer of choice. |
| Storage Instruction | Store at -20°C for long term storage. Avoid freeze-thaw cycles.  |

### TARGET INFORMATION

|             |  |
|-------------|--|
| Gene ID     | <a href="#">227288</a>   |
| Gene Symbol | <a href="#">Cxcr1</a>  |
| UniProt ID  | <a href="#">CXCR1_MOUSE</a>  |
| Specificity | This blocking peptide is recommended for use in combination with CXCR1 antibody, STJ500647 |

### ADDITIONAL INFORMATION

|      |   |
|------|---|
| Note | STRICTLY FOR FURTHER SCIENTIFIC RESEARCH USE ONLY (RUO). MUST NOT TO BE USED IN DIAGNOSTIC OR THERAPEUTIC APPLICATIONS. |
|------|---|