

C6orf150 Positive Control for STJ500507 peptide (STJ503902)

GENERAL INFORMATION

Product Type	Biomolecules
Applications	WB

PRODUCT PROPERTIES

Dilution Range	WB: 1:500
Formulation	Provided as 100 uL ready-to-use, in SDS-PAGE sample buffer (Laemelli's buffer) containing Tris, pH 6.8, 1 % SDS, Glycerol and Bromophenolblue blue as tracking dye. The sample is reduced by adding 2% beta mercaptoethanol. The protein concentration is
Storage Instruction	Store at -20°C for long term storage. Avoid freeze-thaw cycles.

TARGET INFORMATION

Gene ID	115004
Gene Symbol	CGAS
UniProt ID	CGAS_HUMAN
Specificity	This is positive control is recommended for use in combination with C6orf150 antibody STJ500507.

ADDITIONAL INFORMATION

Note	STRICTLY FOR FURTHER SCIENTIFIC RESEARCH USE ONLY (RUO). MUST NOT TO BE USED IN DIAGNOSTIC OR THERAPEUTIC APPLICATIONS.
------	---