

Dopamine D5 Receptor Blocking Peptide for STJ500723 peptide (STJ503945)

GENERAL INFORMATION	
Product Type	Biomolecules
Applications	Immunodepletion/Immunocompetition
PRODUCT PROPERTIES	
Formulation	Liquid form at 2.5mg/ml concentration in PBS. Up to 5% DMSO can be added. Orders with >1mg can be supplied in lyophilized powder form, or in buffer of choice.
Storage Instruction	Store at -20°C for long term storage. Avoid freeze-thaw cycles.
TARGET INFORMATION	
Gene ID	13492
Gene Symbol	Drd5
UniProt ID	DRD5_MOUSE
Immunogen Region	428-478
Specificity	This blocking peptide is recommended for use in combination with Dopamine D5 Receptor antibody, STJ500723
ADDITIONAL INFORMATION	
Note	STRICTLY FOR FURTHER SCIENTIFIC RESEARCH USE ONLY (RUO). MUST NOT TO BE USED IN DIAGNOSTIC OR THERAPEUTIC APPLICATIONS.