

Ephrin Receptor A4 Blocking Peptide for STJ500916 peptide (STJ504012)

GENERAL INFORMATION

Product Type	Biomolecules
Applications	Immunodepletion/Immunocompetition

PRODUCT PROPERTIES

Formulation	Liquid form at 2.5mg/ml concentration in PBS. Up to 5% DMSO can be added. Orders with >1mg can be supplied in lyophilized powder form, or in buffer of choice.
Storage Instruction	Store at -20°C for long term storage. Avoid freeze-thaw cycles.

TARGET INFORMATION

Gene ID	13838
Gene Symbol	Epha4
UniProt ID	EPHA4_MOUSE
Immunogen Region	550-650
Specificity	This blocking peptide is recommended for use in combination with Ephrin Receptor A4 antibody, STJ500916

ADDITIONAL INFORMATION

Note	STRICTLY FOR FURTHER SCIENTIFIC RESEARCH USE ONLY (RUO). MUST NOT TO BE USED IN DIAGNOSTIC OR THERAPEUTIC APPLICATIONS.
------	---