

Ephrin Receptor B3 Blocking Peptide for STJ500938 peptide (STJ504019)

GENERAL INFORMATION

Product Type	Biomolecules
Applications	Immunodepletion/Immunocompetition

PRODUCT PROPERTIES

Formulation	Liquid form at 2.5mg/ml concentration in PBS. Up to 5% DMSO can be added. Orders with >1mg can be supplied in lyophilized powder form, or in buffer of choice.
Storage Instruction	Store at -20°C for long term storage. Avoid freeze-thaw cycles.

TARGET INFORMATION

Gene ID	13845
Gene Symbol	Ephb3
UniProt ID	EPHB3_MOUSE
Immunogen Region	100-200
Specificity	This blocking peptide is recommended for use in combination with Ephrin Receptor B3 antibody, STJ500938

ADDITIONAL INFORMATION

Note	STRICTLY FOR FURTHER SCIENTIFIC RESEARCH USE ONLY (RUO). MUST NOT TO BE USED IN DIAGNOSTIC OR THERAPEUTIC APPLICATIONS.
------	---