

KCNJ5 Blocking Peptide for STJ501172 peptide (STJ504132)

GENERAL INFORMATION	
Product Type	Biomolecules
Applications	Immunodepletion/Immunocompetition
PRODUCT PROPERTIES	
Formulation	Liquid form at 2.5mg/ml concentration in PBS. Up to 5% DMSO can be added. Orders with >1mg can be supplied in lyophilized powder form, or in buffer of choice.
Storage Instruction	Store at -20°C for long term storage. Avoid freeze-thaw cycles.
TARGET INFORMATION	
Gene ID	3762
Gene Symbol	KCNJ5
UniProt ID	KCNJ5_HUMAN
Immunogen Region	350-450
Specificity	This blocking peptide is recommended for use in combination with KCNJ5 antibody, STJ501172
ADDITIONAL INFORMATION	
Note	STRICTLY FOR FURTHER SCIENTIFIC RESEARCH USE ONLY (RUO). MUST NOT TO BE USED IN DIAGNOSTIC OR THERAPEUTIC APPLICATIONS.