

## GRIK1 Blocking Peptide for STJ501192 peptide (STJ504136)

### GENERAL INFORMATION

|              |                                   |
|--------------|-----------------------------------|
| Product Type | Biomolecules                      |
| Applications | Immunodepletion/Immunocompetition |

### PRODUCT PROPERTIES

|                     |                                                                                                                                                                |
|---------------------|----------------------------------------------------------------------------------------------------------------------------------------------------------------|
| Formulation         | Liquid form at 2.5mg/ml concentration in PBS. Up to 5% DMSO can be added. Orders with >1mg can be supplied in lyophilized powder form, or in buffer of choice. |
| Storage Instruction | Store at -20°C for long term storage. Avoid freeze-thaw cycles.                                                                                                |

### TARGET INFORMATION

|             |                                                                                            |
|-------------|--------------------------------------------------------------------------------------------|
| Specificity | This blocking peptide is recommended for use in combination with GRIK1 antibody, STJ501192 |
|-------------|--------------------------------------------------------------------------------------------|

### ADDITIONAL INFORMATION

|      |                                                                                                                         |
|------|-------------------------------------------------------------------------------------------------------------------------|
| Note | STRICTLY FOR FURTHER SCIENTIFIC RESEARCH USE ONLY (RUO). MUST NOT TO BE USED IN DIAGNOSTIC OR THERAPEUTIC APPLICATIONS. |
|------|-------------------------------------------------------------------------------------------------------------------------|