

GPCR GPR3 Blocking Peptide for STJ501255 peptide (STJ504156)

GENERAL INFORMATION	
Product Type	Biomolecules
Applications	Immunodepletion/Immunocompetition
PRODUCT PROPERTIES	
Formulation	Liquid form at 2.5mg/ml concentration in PBS. Up to 5% DMSO can be added. Orders with >1mg can be supplied in lyophilized powder form, or in buffer of choice.
Storage Instruction	Store at -20°C for long term storage. Avoid freeze-thaw cycles.
TARGET INFORMATION	
Gene ID	266769
Gene Symbol	Gpr3
UniProt ID	GPR3_RAT
Immunogen Region	5-50
Specificity	This blocking peptide is recommended for use in combination with GPCR GPR3 antibody, STJ501255
ADDITIONAL INFORMATION	
Note	STRICTLY FOR FURTHER SCIENTIFIC RESEARCH USE ONLY (RUO). MUST NOT TO BE USED IN DIAGNOSTIC OR THERAPEUTIC APPLICATIONS.