

## HIF-1-Alpha isoform 1 + 2 Blocking Peptide for STJ501338 peptide (STJ504234)

### GENERAL INFORMATION

Product Type	Biomolecules
Applications	Immunodepletion/Immunocompetition

### PRODUCT PROPERTIES

Formulation	Liquid form at 2.5mg/ml concentration in PBS. Up to 5% DMSO can be added. Orders with >1mg can be supplied in lyophilized powder form, or in buffer of choice.
Storage Instruction	Store at -20°C for long term storage. Avoid freeze-thaw cycles.

### TARGET INFORMATION

Gene ID	<a href="#">3091</a>
Gene Symbol	<a href="#">HIF1A</a>
UniProt ID	<a href="#">HIF1A_HUMAN</a>
Immunogen Region	150-250
Specificity	This blocking peptide is recommended for use in combination with HIF-1-Alpha isoform 1 + 2 antibody, STJ501338

### ADDITIONAL INFORMATION

Note	STRICTLY FOR FURTHER SCIENTIFIC RESEARCH USE ONLY (RUO). MUST NOT TO BE USED IN DIAGNOSTIC OR THERAPEUTIC APPLICATIONS.
------	---