

IL32 Positive Control for STJ501448 peptide (STJ504317)

GENERAL INFORMATION	
Product Type	Biomolecules
Applications	WB
PRODUCT PROPERTIES	
Dilution Range	WB: 1:500
Formulation	Provided as 100 uL ready-to-use, in SDS-PAGE sample buffer (Laemelli's buffer) containing Tris, pH 6.8, 1 % SDS, Glycerol and Bromophenolblue blue as tracking dye. The sample is reduced by adding 2% beta mercaptoethanol. The protein concentration is
Storage Instruction	Store at -20°C for long term storage. Avoid freeze-thaw cycles.
TARGET INFORMATION	
Gene ID	9235
Gene Symbol	IL32
UniProt ID	IL32_HUMAN
Specificity	This is positive control is recommended for use in combination with IL32 antibody STJ501448.
ADDITIONAL INFORMATION	
Note	STRICTLY FOR FURTHER SCIENTIFIC RESEARCH USE ONLY (RUO). MUST NOT TO BE USED IN DIAGNOSTIC OR THERAPEUTIC APPLICATIONS.