

MUC1 Blocking Peptide for STJ501818 peptide (STJ504390)

GENERAL INFORMATION

Product Type	Biomolecules
Applications	Immunodepletion/Immunocompetition

PRODUCT PROPERTIES

Dilution Range	WB: IHC: 1:200 ICC: 1:200 IF: 1:200 CFM: 1:200
Formulation	Liquid form at 2.5mg/ml concentration in PBS. Up to 5% DMSO can be added. Orders with >1mg can be supplied in lyophilized powder form, or in buffer of choice.
Storage Instruction	Store at -20°C for long term storage. Avoid freeze-thaw cycles.

TARGET INFORMATION

Gene ID	4582
Gene Symbol	MUC1
UniProt ID	MUC1_HUMAN
Specificity	This blocking peptide is recommended for use in combination with MUC1 antibody, STJ501818

ADDITIONAL INFORMATION

Note	STRICTLY FOR FURTHER SCIENTIFIC RESEARCH USE ONLY (RUO). MUST NOT TO BE USED IN DIAGNOSTIC OR THERAPEUTIC APPLICATIONS.
------	---