

## MyoD1 Blocking Peptide peptide (STJ504397)

### GENERAL INFORMATION

|              |                                   |
|--------------|-----------------------------------|
| Product Type | Biomolecules                      |
| Applications | Immunodepletion/Immunocompetition |

### PRODUCT PROPERTIES

|                     |                                                                                                                                                                |
|---------------------|----------------------------------------------------------------------------------------------------------------------------------------------------------------|
| Formulation         | Liquid form at 2.5mg/ml concentration in PBS. Up to 5% DMSO can be added. Orders with >1mg can be supplied in lyophilized powder form, or in buffer of choice. |
| Storage Instruction | Store at -20°C for long term storage. Avoid freeze-thaw cycles.                                                                                                |

### TARGET INFORMATION

|                  |                             |
|------------------|-----------------------------|
| Gene ID          | <a href="#">4654</a>        |
| Gene Symbol      | <a href="#">MYOD1</a>       |
| UniProt ID       | <a href="#">MYOD1_HUMAN</a> |
| Immunogen Region | 280-330                     |

### ADDITIONAL INFORMATION

|      |                                                                                                                         |
|------|-------------------------------------------------------------------------------------------------------------------------|
| Note | STRICTLY FOR FURTHER SCIENTIFIC RESEARCH USE ONLY (RUO). MUST NOT TO BE USED IN DIAGNOSTIC OR THERAPEUTIC APPLICATIONS. |
|------|-------------------------------------------------------------------------------------------------------------------------|