

NPY1R Blocking Peptide for STJ501961 peptide (STJ504426)

GENERAL INFORMATION

| | |
|--------------|-----------------------------------|
| Product Type | Biomolecules |
| Applications | Immunodepletion/Immunocompetition |

PRODUCT PROPERTIES

| | |
|---------------------|--|
| Formulation | Liquid form at 2.5mg/ml concentration in PBS. Up to 5% DMSO can be added. Orders with >1mg can be supplied in lyophilized powder form, or in buffer of choice. |
| Storage Instruction | Store at -20°C for long term storage. Avoid freeze-thaw cycles. |

TARGET INFORMATION

| | |
|------------------|--|
| Immunogen Region | 360-410 |
| Specificity | This blocking peptide is recommended for use in combination with NPY1R antibody, STJ501961 |

ADDITIONAL INFORMATION

| | |
|------|---|
| Note | STRICTLY FOR FURTHER SCIENTIFIC RESEARCH USE ONLY (RUO). MUST NOT TO BE USED IN DIAGNOSTIC OR THERAPEUTIC APPLICATIONS. |
|------|---|