

alpha Actinin 4 Blocking Peptide peptide (Phospho) (STJ504522)

GENERAL INFORMATION

Product Type	Biomolecules
Applications	Immunodepletion/Immunocompetition

PRODUCT PROPERTIES

Storage Instruction	Store at -20°C for long term storage. Avoid freeze-thaw cycles.
---------------------	---

TARGET INFORMATION

Gene ID	81
Gene Symbol	ACTN4
UniProt ID	ACTN4_HUMAN
Immunogen Region	60-110

ADDITIONAL INFORMATION

Note	STRICTLY FOR FURTHER SCIENTIFIC RESEARCH USE ONLY (RUO). MUST NOT TO BE USED IN DIAGNOSTIC OR THERAPEUTIC APPLICATIONS.
------	---