

Chimaerin beta 2 Positive Control peptide (Phospho) (STJ504584)

GENERAL INFORMATION

Product Type	Biomolecules
Applications	WB

PRODUCT PROPERTIES

Dilution Range	WB: 1:500
Storage Instruction	Store at -20°C for long term storage. Avoid freeze-thaw cycles.

TARGET INFORMATION

Gene ID	1123
Gene Symbol	CHN1
UniProt ID	CHIN_HUMAN

ADDITIONAL INFORMATION

Note	STRICTLY FOR FURTHER SCIENTIFIC RESEARCH USE ONLY (RUO). MUST NOT TO BE USED IN DIAGNOSTIC OR THERAPEUTIC APPLICATIONS.
------	---