

PDE3A Blocking Peptide for STJ502266 peptide (STJ504612)

GENERAL INFORMATION

Product Type Biomolecules

PRODUCT PROPERTIES

Formulation Liquid form at 2.5mg/ml concentration in PBS. Up to 5% DMSO can be added. Orders with >1mg can be supplied in lyophilized powder form, or in buffer of choice.

Storage Instruction Store at -20°C for long term storage. Avoid freeze-thaw cycles.

TARGET INFORMATION

Gene ID [5139](#)

Gene Symbol [PDE3A](#)

UniProt ID [PDE3A_HUMAN](#)

Specificity This blocking peptide is recommended for use in combination with PDE3A antibody, STJ502266

ADDITIONAL INFORMATION

Note STRICTLY FOR FURTHER SCIENTIFIC RESEARCH USE ONLY (RUO). MUST NOT TO BE USED IN DIAGNOSTIC OR THERAPEUTIC APPLICATIONS.